



Adenovirus antigen (ADV Ag) Rapid Tests (DST-H008S)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Intended Use	The CD ADV Ag Rapid Test is a lateral flow immunochromatographic assay for the qualitative detection of Adenovirus antigen in human stool. It is intended to be used as a screening test and as an aid in the diagnosis of infection with Adenovirus.
Reagents And Materials Provided	For strip test: 20 X Test devices. Each test strp is packed in a foil pouch with a package of desiccant. 20 X Assay buffer (Stool buffer tubes). 1 X Instruction for use.
Storage	Store at 2-30°C. Expiry period 24 months.
Performance Characteristics	1. State-of-arts workmanships and severe lab testing processes; Product CE approved and quality system ISO13485 certified; 2. Easy assay procedure and rapid results reading; High sensitivity and specificity; 3. Full-aspects OEM service provided. Design according to customers unique requirements. Help improving customer's brand value.
Precision	98.5%
Sensitivity	95.6%