



Clenbuterol Rapid Test (Feed) (DTS419)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Intended Use	CD Clenbuterol Rapid Test is a competitive immunoassay for the semi-quantitative detection of the presence of Clenbuterol residue in feed.
Storage	The kit can be stored at room temperature (2-30°C). The test kit is stable through the expiration date (18 months) marked on the foil pouch. DO NOT FREEZE. Do not store the test kit in direct sunlight.
Sensitivity	10 ppb
Specificity	The results are negative when the test is applied to detect 100 ppm of Ractopamine, Trenbolone-acetate, Chloramphenicol, Aminoglycosides, Beta-Lactams, Macrolides, Quinolones, and Sulfonamides.