



Ractopamine Rapid Test (Tissue) (DTS422)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Intended Use	CD Ractopamine Rapid Test is a competitive immunoassay for the semi-quantitative detection of the presence of Ractopamine residue in animal's tissue.
Storage	The kit can be stored at room temperature (2-30°C). The test kit is stable through the expiration date (18 months) marked on the foil pouch. DO NOT FREEZE. Do not store the test kit in direct sunlight.
Sensitivity	5 ppb
Specificity	The results are negative when the test is applied to detecting 100 ppm of Clenbuterol, Salbutamol, Trenbolone-acetate, Chloramphenicol, Aminoglycosides, Beta-Lactams, Macrolides, Quinolones, and Tetracyclines.