



Aflatoxin B1 Precise Test (Fermented Wine) (DTS433)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Intended Use	CD Aflatoxin B1 Precise Test is a competitive immunoassay for the semi-quantitative detection of the presence of Aflatoxin B1 in fermented wine.
Storage	The kit can be stored at room temperature (2-30°C). The test kit is stable through the expiration date (18 months) marked on the foil pouch. DO NOT FREEZE. Do not store the test kit in direct sunlight.
Sensitivity	5 ppb