



## Deoxynivalenol (Vomitoxin) Test (DTS436)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Intended Use</b>	CD Deoxynivalenol Rapid Test is a competitive immunoassay for the semi-quantitative detection of the presence of Deoxynivalenol in grain or feed.
<b>Storage</b>	The kit can be stored at room temperature (2-30°C). The test kit is stable through the expiration date (18 months) marked on the foil pouch. DO NOT FREEZE. Do not store the test kit in direct sunlight.