



Salbutamol Rapid Test (Tissue) (DTS443)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Intended Use	CD Salbutamol Rapid Test is a useful tool offering a rapid and accurate testing in field screening, exceeding with its convenience. It provides a semi-quantitative method to detect the Salbutamol above 5 ppb in meat. If you want a quantitative result, please adopt other method such as ELISA in practice.
Storage	The kit can be stored at room temperature (2-30°C). The test kit is stable through the expiration date (18 months) marked on the foil pouch. DO NOT FREEZE. Do not store the test kit in direct sunlight.
Sensitivity	5 ppb