



Sulfadimethoxine Rapid Test (Honey) (DTS445)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Intended Use	CD Sulfadimethoxine Rapid Test is a competitive immunoassay for the semi-quantitative detection of the presence of Sulfadimethoxine residue in honey sample.
	Cut-off : 10 ppb
	Assay Time : 10-15 min
Storage	<p>The kit can be stored at room temperature (2-30°C). The test kit is stable through the expiration date (18 months) marked on the foil pouch.</p> <p>DO NOT FREEZE</p> <p>. Do not store the test kit in direct sunlight.</p>
Sensitivity	10 ppb
Specificity	The results are negative when the test card is applied to detecting 100 ppm (µg/mL) of Chloramphenicol, Aminoglycosides, Beta-lactams, Fluoro quinolones, and Macrolides.