



## Ceftiofur Rapid Test (Milk) (DTS448)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

#### Intended Use

CD Ceftiofur Rapid Test is a competitive immunoassay for the semi-quantitative detection of the presence of Ceftiofur residue in milk.

Cut-off  
: 50 ppb

Assay Time  
: 5-10 min

Sample  
: Defatted milk

---

#### Storage

The kit can be stored at room temperature (2-30°C). The test kit is stable through the expiration date (18 months) marked on the foil pouch.

DO NOT FREEZE

. Do not store the test kit in direct sunlight.

---

<b>Sensitivity</b>	50 ppb
--------------------	--------

---

<b>Specificity</b>	The results are negative when the test card is applied to detect 100 ppm of Chloramphenicol, Macrolides, Sulfonamides, Cephapirin, Aminoglycosides, Tetracycline and Penicillin.
--------------------	--

---