



Mouse Anti-Metalaxyl monoclonal antibody, clone MA (HMABPY0003)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Metalaxyl with carrier protein.
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	MA
Purification	Purified from mouse ascites.
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	100 µg, 1 mg
Buffer	PBS
Preservative	None
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Introduction

Metalaxyl is an acylalanine fungicide with systemic function. Its chemical name is methyl N-(methoxyacetyl)-N-(2,6-xylyl)-DL-alaninate. It can be used to control Pythium in a number of vegetable crops, and Phytophthora in peas. Metalaxyl-M or Ridomil Gold are trade names for the optically pure (-) / D / R active stereoisomer, which is also known as Mefenoxam.

Keywords

Metalaxyl;METALAXYL;57837-19-1;Ridomil;Allegiance;Metanaxin;