



# RHA™ anti-Ractopamine monoclonal antibody, clone RCT (HMABPY001)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

### Antigen Description

Ractopamine is a phenolethanolamine beta-adrenoceptor agonist that is used for the improvement of weight gain, carcass leanness and feed efficiency in pigs. It has been known for years that ractopamine enhances growth in meat animals by repartitioning nutrients. This means that the nutrients that the animal consumes in its feed are shifted away from being stored as fat and are instead used in lean tissue production. This shift occurs as the animal reaches the top of its growth curve and lean tissue production drops off and fat deposition increases. Fat levels also appear to decrease, due to an increase in lipolysis (fat breakdown). Using the ractopamine monoclonal antibody, it is possible to detect ractopamine and other agonists in urine, muscle and liver both rapidly and with accuracy.

<b>Immunogen</b>	Ractopamine [KLH]
<b>Isotype</b>	IgG
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	N/A
<b>Clone</b>	RCT
<b>Purification</b>	Purified from mouse ascites.
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, LFIA
<b>Size</b>	100 µg, 1 mg
<b>Buffer</b>	0.01M pH7.4 PBS
<b>Preservative</b>	None

**Storage**

Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

---