



RHA™ anti-Clenbuterol residues monoclonal antibody, clone β 2-AR (HMABPY004)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Clenbuterol Residues [KLH]
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	β 2-AR
Purification	Purified from mouse ascites.
Conjugate	Unconjugated
Applications	ELISA, LFIA
Size	100 μ g, 1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.