



RHA™ anti-Zearalenone monoclonal antibody, clone ZEN (HMABPY013)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Cross reactivity: Zearalanone: 0.04ppb, β -Zearalanol: 0.08ppb, β -Zearalenol: 0.07ppb, α -Zearalanol: 0.04ppb, α -Zearalenol: 0.06ppb
Immunogen	Zearalenone [KLH]
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	ZEN
Purification	Purified from mouse ascites.
Conjugate	Unconjugated
Applications	ELISA, LFIA
Size	100 μ g, 1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.