



RHA™ anti-Thiamphenicol drug residue monoclonal antibody, clone tap (HMABPY041)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|--------------------|---|
| Immunogen | Thiamphenicol Drug [KLH] |
| Isotype | IgG |
| Source/Host | Mouse |
| Species Reactivity | N/A |
| Clone | tap |
| Purification | Purified from mouse ascites. |
| Conjugate | Unconjugated |
| Applications | ELISA, LFIA |
| Size | 100 µg, 1 mg |
| Buffer | 0.01M pH7.4 PBS |
| Preservative | None |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |