



RHA™ anti-Avermectin monoclonal antibody, clone AVM (HMABPY051)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Cross-Reactivity: Enamectin 0.5ppb; Ivermectin 1.5ppb; Doramectin 20ppb
Immunogen	Avermectin [KLH]
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	AVM
Purification	Purified from mouse ascites.
Conjugate	Unconjugated
Applications	ELISA, LFIA
Size	100 µg, 1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.