



## RHA™ anti-Albendazole monoclonal antibody, clone abz (HMABPY057)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Specificity</b>	Cross-Reactivity: ABZSO <sub>2</sub> NH <sub>2</sub> 10ppb; ABZSO <sub>2</sub> 10ppb; ABZSO 80ppb
<b>Immunogen</b>	Albendazole [KLH]
<b>Isotype</b>	IgG
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	N/A
<b>Clone</b>	abz
<b>Purification</b>	Purified from mouse ascites.
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, LFIA
<b>Size</b>	100 µg, 1 mg
<b>Buffer</b>	0.01M pH7.4 PBS
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.