



# RHA™ anti-Fenpropathrin monoclonal antibody, clone FPP (HMABPY092)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	Fenpropathrin (brand names Danitol, Meothrin), or fenopropathrin, is a widely used pyrethroid insecticide in agriculture and household. It has thus been implicated as an environmental risk factor for Parkinson's disease.
<b>Immunogen</b>	Fenpropathrin [KLH]
<b>Isotype</b>	IgG
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	N/A
<b>Clone</b>	FPP
<b>Purification</b>	Purified from mouse ascites.
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, LFIA
<b>Size</b>	100 µg, 1 mg
<b>Buffer</b>	0.01M pH7.4 PBS
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.