



RHA™ anti-Melamine monoclonal antibody, clone MEL (HMABPY096)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Cross-Reactivity: Cyromazine 2ppb
Immunogen	Melamine [KLH]
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	MEL
Purification	Purified from mouse ascites.
Conjugate	Unconjugated
Applications	ELISA, LFIA
Size	100 µg, 1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.