



RHA™ anti-Dibutyl phthalate monoclonal antibody, clone dbp (HMABPY111)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description Dibutyl phthalate (DBP) is a commonly used plasticizer. It is also used as an additive to adhesives or printing inks. It is soluble in various organic solvents, e.g. in alcohol, ether and benzene. DBP is also used as an ectoparasiticide.

Immunogen Dibutyl Phthalate [KLH]

Isotype IgG

Source/Host Mouse

Species Reactivity N/A

Clone dbp

Purification Purified from mouse ascites.

Conjugate Unconjugated

Applications ELISA, LFIA

Size 100 µg, 1 mg

Buffer 0.01M pH7.4 PBS

Preservative None

Storage Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.