



Mouse Anti-Peanut allergen Ara H2 monoclonal antibody, clone Arah2D (HMABPY117D)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Peanut allergen Ara H2
Target	Peanut allergen Ara H2
Immunogen	Peanut allergen Ara H2
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Arachis hypogaea
Clone	Arah2D
Purification	Purified from mouse ascites.
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	100 µg, 1 mg
Buffer	PBS
Preservative	None
Storage	Long time storage is recommended at -20°C.

Ship

Wet ice

BACKGROUND

Introduction

Ara h 2 is a glycoprotein purified from light roasted peanut flour by multi-step chromatography. Suitable for research studies of food allergy.

Keywords

Ara h 2;Peanut allergen Ara H2
