



One Step Malaria Pf/Pan (Whole blood) Test (IVDTS003)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Intended Use

One Step Malaria Pf/Pan Test is a lateral flow chromatographic immunoassay for the detection of Malaria plasmodium lactate dehydrogenase (pLDH) in vitro assay. It is used for detection of Malaria *P. falciparum* and other Plasmodium (*P. vivax*, *P. malariae*, *P. ovale*) in human blood. This device is intended to be used as a screening test and as an aid in the diagnosis of infection with Plasmodium. Any reactive specimen with One Step Malaria Pf/Pan Test must be confirmed with alternative testing method(s) and clinical findings.

Storage

The test kit can be stored at 4 to 30°C in the sealed pouch to the date of expiration. The test kit should be kept away from direct sunlight, moisture and heat.