



# Human HMGB1 ELISA Matched Antibody Pair (ABPR-0429)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

### Intended Use

This antibody pair set comes with matched antibody pair to detect and quantify protein level of human HMGB1.

### Reagents And Materials Provided

Antibody pair set content:

1. Capture antibody: mouse monoclonal anti-HMGB1, IgG2a Kappa (100 µg)
2. Detection antibody: biotinylated mouse monoclonal anti-HMGB1, IgG1 Kappa (50 µg)

Note: Reagents are sufficient for at least 3-5 x 96 well plates using recommended protocols.

### Reconstitution And Storage

Store reagents of the antibody pair set at -20°C or lower. Please aliquot to avoid repeated freeze thaw cycle. Reagents should be returned to -20°C storage immediately after use.