



# Anti-ACBD3 (Internal) polyclonal antibody (CABT-B361)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Target</b>	ACBD3
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse, Dog, Rat, Zebrafish, Yeast, Guinea pig, Horse, Rabbit, Bovine
<b>Purification</b>	Affinity Purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB
<b>Format</b>	Liquid
<b>Concentration</b>	0.5 mg/ml (Please refer to the vial label for the specific concentration)
<b>Size</b>	100 µl
<b>Buffer</b>	In 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
<b>Storage</b>	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.

## BACKGROUND

**Keywords** ACBD3;acyl-CoA binding domain containing 3;PAP7;GCP60;GOCAP1;GOLPH1;Golgi resident protein GCP60;golgi phosphoprotein 1;PKA (RIalpha)-associated protein;PBR- and PKA-associated protein 7;Acbd3

## GENE INFORMATION

**Synonyms** acyl-CoA binding domain containing 3; Acbd3