



Anti-Commd2 polyclonal antibody (CABT-B486)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Commd2
Immunogen	Synthetic peptide corresponding to a region of Rat
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Dog, Pig, Rat, Zebrafish, Guinea pig, Horse, Rabbit, Bovine
Conjugate	Unconjugated
Applications	WB
Format	Lyophilized powder
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration)
Size	50 µg
Buffer	PBS buffer with 2% sucrose.
Preservative	None
Storage	For longer periods of storage, store at -20°C. Avoid repeat freeze-thaw cycles.

BACKGROUND

Keywords	COMMD2;COMM domain containing 2;COMM domain-containing protein 2;HSPC042;MGC57611
-----------------	---

GENE INFORMATION

Synonyms	COMM domain containing 2
-----------------	--------------------------