



Anti-ABF1 (aa 201-300) polyclonal antibody (CABT-BL160)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Fusion protein: RGLDDESGPTHGNDSGNHRHNEEDDVHTQMTKNYSDEVVND EDINVAIANAVANVDSQSNNKHDGKDDDATNNNDGQDNNT NNDHNNNSNINNNNVGSHGI, corresponding to amino acids 201-300 of <i>S. cerevisiae</i> ABF1.
Isotype	IgG
Source/Host	Mouse
Species Reactivity	<i>S. cerevisiae</i>
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nuclear
Format	Liquid
Buffer	50% Glycerol
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	853748
Protein Refseq	NP_012810
UniProt ID	P14164