



# Anti-GPR137C (internal region) polyclonal antibody (CABT-BL1659)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Synthetic peptide derived from an internal region of human GPR137C.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, ELISA, IHC-P, ICC/IF
<b>Cellular Localization</b>	Membrane.
<b>Format</b>	Liquid
<b>Buffer</b>	50% Glycerol, PBS (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), 150mM Sodium chloride, pH 7.4
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Store at -20°C. Stable for 12 months at -20°C

## GENE INFORMATION

<b>Entrez Gene ID</b>	<a href="#">283554</a>
<b>Protein Refseq</b>	<a href="#">NP_001093122</a>
<b>UniProt ID</b>	<a href="#">B3KW22</a>