



Anti-KLRA8 monoclonal antibody, clone 4E20 (CABT-BL2088)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Mouse
Clone	4E20
Conjugate	Unconjugated
Applications	FC, FuncS, IP
Preservative	None
Storage	2-8°C (up to 1 month) below -15°C (long term)

GENE INFORMATION

Entrez Gene ID	16639
Protein Refseq	NP_034780.1
UniProt ID	Q60682