



# Anti-KLRA8 monoclonal antibody, clone 4E20 [Alexa Fluor® 647] (CABT-BL2090)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Isotype</b>	IgG1
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	Mouse
<b>Clone</b>	4E20
<b>Conjugate</b>	Alexa Fluor® 647
<b>Applications</b>	FC
<b>Preservative</b>	None
<b>Storage</b>	2-8°C (up to 1 month) below -15°C (long term)

## GENE INFORMATION

<b>Entrez Gene ID</b>	<a href="#">16639</a>
<b>Protein Refseq</b>	<a href="#">NP_034780.1</a>
<b>UniProt ID</b>	<a href="#">Q60682</a>