



Anti-PAIN (cytoplasmic domain) polyclonal antibody (CABT-BL2847)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide from 4th cytoplasmic loop of drosophila Pain conjugated to an immunogenic carrier protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Fruit fly
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB, IHC-P, IHC-Fr
Reconstitution	Reconstitute in 100 µl of sterile water. Centrifuge to remove any insoluble material.
Cellular Localization	Membrane; Multi-pass membrane protein
Format	Lyophilised
Buffer	Whole serum
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C (add glycerol to a final volume of 40% for extra stability). Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	37985
Protein Refseq	NP_611979
UniProt ID	Q9W0Y6