



Anti-PUT (aa 128-147) polyclonal antibody (CABT-BL3066)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide: KSTTEATTQVPKEKTQDGSN, corresponding to C terminal amino acids 128-147 of <i>Drosophila melanogaster</i> Punt.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>Drosophila</i>
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IP, WB, ICC
Format	Liquid
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	41772
Protein Refseq	NP_731926
UniProt ID	Q24468