



Anti-ARP5 (full length) polyclonal antibody (CABT-BL341)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant full length protein Arp 5.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>S. cerevisiae</i>
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nuclear
Format	Liquid
Buffer	1% BSA, PBS, pH 7.4
Preservative	0.02% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	855665
Protein Refseq	NP_014339
UniProt ID	P53946