



## Anti-TAS2R41 polyclonal antibody (CABT-BL3524)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Immunogen</b>	A synthetic peptide as a part of human TAS2R41 conjugated to immunogenic carrier protein
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Whole antiserum
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, IHC-P, IHC-Fr
<b>Cellular Localization</b>	Membrane; Multi-pass membrane protein.
<b>Format</b>	Liquid
<b>Buffer</b>	PBS, Whole serum
<b>Preservative</b>	None
<b>Storage</b>	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C (add glycerol to a final volume of 40% for extra stability). Avoid repeated freeze / thaw cycles.

### GENE INFORMATION

<b>Entrez Gene ID</b>	<a href="#">259287</a>
<b>Protein Refseq</b>	<a href="#">NP_795364</a>
<b>UniProt ID</b>	<a href="#">P59536</a>