



Anti-BEND2 (full length) polyclonal antibody (CABT-BL4216)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse Polyclonal antibody to Human BEND2.
Immunogen	Full length human CXORF20 protein.
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Purification	Protein G purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	1×PBS, pH 7.2
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze/thaw cycles.

BACKGROUND

Introduction	This gene encodes a protein which has two BEN domains in the C-terminus. These domains are found in proteins which participate in protein and DNA interactions which occur during chromatin restructuring or transcription. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2011]
---------------------	--

GENE INFORMATION

Entrez Gene ID	139105
Protein Refseq	NP_001171696
UniProt ID	Q8NDZ0
Chromosome Location	Xp22.22