



Anti-PACRGL polyclonal antibody (CABT-BL5341)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-PACRGL polyclonal Antibody
Target	PACRGL
Immunogen	Recombinant protein corresponding to amino acids of human PACRGL.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Sequence Similarities	ARKLHPRPSDKLNPKTINPFGEQSRVPSAFAAIYSKGGIPCRLVHGSKHRLQWECPPESLSFD PLLITLAEGLR ETKHPYTFVSKEGFRELLLVKGAPE
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Entrez Gene ID	133015
Protein Refseq	NP_001124199
UniProt ID	Q8N7B6
Chromosome Location	4p15.31
Function	binding;