



Anti-PAMR1 (internal region) polyclonal antibody (CABT-BL5355)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit polyclonal antibody to Human PAMR1.
Immunogen	Synthetic peptide derived from the internal domain of Human PAMR1 protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Whole antiserum
Conjugate	Unconjugated
Cellular Localization	Secreted
Format	Liquid
Size	200 µl
Buffer	Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze/thaw cycles.

GENE INFORMATION

Entrez Gene ID	25891
Protein Refseq	NP_056245
UniProt ID	Q6UXH9
Chromosome Location	11p13
Function	catalytic activity; serine-type endopeptidase activity;