



Anti-PCBD2 (full length) polyclonal antibody (CABT-BL5359)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse polyclonal antibody to Human PCBD2.
Immunogen	PCBD2 full-length human protein.
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Purification	Whole antiserum
Conjugate	Unconjugated
Sequence Similarities	Belongs to the pterin-4-alpha-carbinolamine dehydratase family.
Format	Liquid
Size	50 µl
Buffer	Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze/thaw cycles.

GENE INFORMATION

Entrez Gene ID	84105
Protein Refseq	NP_115527
UniProt ID	Q9H0N5
Chromosome Location	5q31.1
Function	4-alpha-hydroxytetrahydrobiopterin dehydratase activity; lyase activity; phenylalanine 4-monooxygenase activity; protein binding;