



Anti-PLEKHN1 polyclonal antibody (CABT-BL5420)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-PLEKHN1 polyclonal Antibody
Target	PLEKHN1
Immunogen	Recombinant protein corresponding to amino acids of human PLEKHN1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Sequence Similarities	IFSEELDGLCFKGELPLRAVHINLEEKEKQIRSFLIEGPLINTIRVVCASYEDYGHWLLCLRAVTHR EGAPPLPG AESFPGSQV
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Entrez Gene ID	84069
Protein Refseq	NP_115505.2
UniProt ID	Q494U1
Chromosome Location	1p36.33
Function	phospholipid binding;