



# Anti-TAS2R60 polyclonal antibody (CABT-BL5818)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Rabbit polyclonal antibody to Human TAS2R60.
<b>Immunogen</b>	A synthetic peptide as a part of human T2R60 conjugated to immunogenic carrier protein.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Whole antiserum
<b>Conjugate</b>	Unconjugated
<b>Cellular Localization</b>	Cell Membrane; Multi-pass membrane protein
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	PBS, Whole serum
<b>Preservative</b>	None
<b>Storage</b>	Store at 4°C short term (1-2 weeks). Aliquot and store at -20°C (add glycerol to a final volume of 40% for extra stability). Avoid repeated freeze/thaw cycles.

## GENE INFORMATION

<b>Entrez Gene ID</b>	<a href="#">338398</a>
<b>Protein Refseq</b>	<a href="#">NP_803186</a>
<b>UniProt ID</b>	<a href="#">P59551</a>
<b>Chromosome Location</b>	7
<b>Pathway</b>	Taste transduction, organism-specific biosystem; Taste transduction, conserved biosystem;
<b>Function</b>	G-protein coupled receptor activity; receptor activity; signal transducer activity;