



Anti-ACLSV polyclonal antibody (CABT-BL6125)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Antiserum is developed in rabbit using purified virus as immunogen.
Source/Host	Rabbit
Purification	caprylic acid purified
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	0.2 ml/1 ml
Buffer	In 0.1 M phosphate buffered saline, pH 7.4.
Storage	The antibody should be stored at 2-8°C. For storage longer than one year, the solution may be frozen at -20°C. Repeated freezing and thawing is not recommended. The solution may be frozen in aliquots if necessary. The antibody has a shelf-life of 2 years after date of purchase.

BACKGROUND

Introduction	A virus with flexuous filamentous particles 680-780 nm long and 12 nm in width containing linear positive-sense, ssRNA. Occurs in woody plants of the family Rosaceae including apple, peach, pear, plum, cherry, apricot and prune. Its natural mode of spread is unknown. Transmissible by mechanical inoculation to herbaceous plants.
Keywords	Apple chlorotic leafspot virus; ACLSV; Pear ring pattern mosaic virus; apple latent virus type 1; plum pseudopox virus; quince stunt virus