



Anti-CNFV polyclonal antibody (CABT-BL6134)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Antiserum is developed in rabbit using purified virus as immunogen.
Isotype	IgG
Source/Host	Rabbit
Purification	caprylic acid purified
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	0.2 ml/1 ml
Buffer	In 0.1 M phosphate buffered saline, pH 7.4.
Preservative	None
Storage	The antibody should be stored at 2-8°C. For storage longer than one year, the solution may be frozen at -20°C. Repeated freezing and thawing is not recommended. The solution may be frozen in aliquots if necessary. The antibody has a shelf-life of 2 years after date of purchase.

BACKGROUND

Introduction	A virus with very flexuous elongated particles c. 1400-1500 nm long. Transmitted by aphids in a semi-persistent manner but with difficulty by inoculation of sap. It seems to infect only species of Caryophyllaceae. Found in Japan.
Keywords	Carnation necrotic fleck virus; CNFV; carnation yellow fleck virus; carnation streak virus