



Anti-SLRSV polyclonal antibody (CABT-BL6173)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Antiserum is developed in rabbit using purified virus as immunogen.
Source/Host	Rabbit
Purification	caprylic acid purified
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	0.2 ml/1 ml
Buffer	In 0.1 M phosphate buffered saline, pH 7.4.
Storage	The antibody should be stored at 2-8°C. For storage longer than one year, the solution may be frozen at -20°C. Repeated freezing and thawing is not recommended. The solution may be frozen in aliquots if necessary. The antibody has a shelf-life of 2 years after date of purchase.

BACKGROUND

Introduction	Strawberry latent ringspot virus (SLRSV) is a plant pathogenic virus in an Unassigned family. Its entry in the ICTV database indicates this virus is a Sadwavirus but current literature indicates that it may be a Cheravirus.
Keywords	Strawberry latent ringspot virus; SLRSV; aesculus line pattern virus; rhubarb virus 5