



Anti-TMGMV polyclonal antibody (CABT-BL6174)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Antiserum is developed in rabbit using purified virus as immunogen.
Source/Host	Rabbit
Purification	caprylic acid purified
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	0.2 ml/1 ml
Buffer	In 0.1 M phosphate buffered saline, pH 7.4.
Storage	The antibody should be stored at 2-8°C. For storage longer than one year, the solution may be frozen at -20°C. Repeated freezing and thawing is not recommended. The solution may be frozen in aliquots if necessary. The antibody has a shelf-life of 12 months.

BACKGROUND

Introduction	A virus with RNA-containing tubular rod-shaped particles c. 308 x 18 nm, easily transmitted by mechanical inoculation and by handling during cultivation. No natural vector known. Widely distributed throughout the world. Of economic importance for tobacco crops.
Keywords	Tobacco mild green mosaic virus; TMGMV; Mild dark-green tobacco mosaic virus; Para-tobacco mosaic virus; Mild strain of tobacco mosaic virus; South Carolina mild mottling strain of tobacco mosaic virus; Strains U2 and U5 of tobacco mosaic virus