



Anti-AICV polyclonal antibody [AP] (CABT-BL6216)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Antiserum is developed in rabbit using purified virus as immunogen.
Source/Host	Rabbit
Purification	caprylic acid purified
Conjugate	AP
Applications	ELISA
Format	Liquid
Size	0.2 ml/1 ml
Buffer	In 0.05 M TRIS buffer, pH 8.0.
Storage	The antibody should be stored at 2-8°C. For storage longer than one year, the solution may be frozen at -20°C. Repeated freezing and thawing is not recommended. The solution may be frozen in aliquots if necessary. The antibody has a shelf-life of 2 years after date of purchase.

BACKGROUND

Introduction	Alstroemeria carlavirus is a virus which can infect Alstroemeria plants. Alstroemeria is an important cut flower. It is mainly cultivated in Netherlands, Columbia, England and America. Alstroemeria is propagated through the division of rhizomes and once the virus infected material is propagated then the virus can transfer from one generation to the next.
Keywords	Alstroemeria carlavirus; AICV