



# Anti-AIFBV polyclonal antibody [AP] (CABT-BL6217)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Antiserum is developed in rabbit using purified virus as immunogen.
<b>Source/Host</b>	Rabbit
<b>Purification</b>	caprylic acid purified
<b>Conjugate</b>	AP
<b>Applications</b>	ELISA
<b>Format</b>	Liquid
<b>Size</b>	0.2 ml/1 ml
<b>Buffer</b>	In 0.05 M TRIS buffer, pH 8.0.
<b>Preservative</b>	None
<b>Storage</b>	The antibody should be stored at 2-8°C. For storage longer than one year, the solution may be frozen at -20°C. Repeated freezing and thawing is not recommended. The solution may be frozen in aliquots if necessary. The antibody has a shelf-life of 2 years after date of purchase.

## BACKGROUND

<b>Introduction</b>	AIFBV is a virus which can infect Alstroemeria plants and causes leaf chlorosis and some flower breakage.
<b>Keywords</b>	Alstroemeria flower-banding virus; AIFBV