



Anti-CVB polyclonal antibody [AP] (CABT-BL6234)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Antiserum is developed in rabbit using purified virus as immunogen.
Source/Host	Rabbit
Purification	caprylic acid purified
Conjugate	AP
Applications	ELISA
Format	Liquid
Size	0.2 ml/1 ml
Buffer	In 0.05 M TRIS buffer, pH 8.0.
Storage	The antibody should be stored at 2-8°C. For storage longer than one year, the solution may be frozen at -20°C. Repeated freezing and thawing is not recommended. The solution may be frozen in aliquots if necessary. The antibody has a shelf-life of 2 years after date of purchase.

BACKGROUND

Introduction	A virus with rather straight rod-shaped particles about 685 nm long. It has a limited host range, is transmitted by aphids in the non-persistent manner and by inoculation of sap. Widespread throughout the world in cultivated varieties of chrysanthemum.
Keywords	Chrysanthemum virus B; CVB; Chrysanthemum mild mosaic virus