



Anti-MiLBVV polyclonal antibody [AP] (CABT-BL6248)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Antiserum is developed in rabbit using purified virus as immunogen.
Source/Host	Rabbit
Purification	caprylic acid purified
Conjugate	AP
Applications	ELISA
Format	Liquid
Size	0.2 ml/1 ml
Buffer	In 0.05 M TRIS buffer, pH 8.0.
Storage	The antibody should be stored at 2-8°C. For storage longer than one year, the solution may be frozen at -20°C. Repeated freezing and thawing is not recommended. The solution may be frozen in aliquots if necessary. The antibody has a shelf-life of 12 months.

BACKGROUND

Introduction	Lettuce big-vein associated virus (LBVaV) is a pathogenic plant virus. The virus is the cause of lettuce big-vein disease.
Keywords	Mirafiori lettuce big vein virus; MiLBVV