



Anti-PFBV polyclonal antibody [AP] (CABT-BL6253)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Antiserum is developed in rabbit using purified virus as immunogen.
Source/Host	Rabbit
Purification	caprylic acid purified
Conjugate	AP
Applications	ELISA
Format	Liquid
Size	0.2 ml/1 ml
Buffer	In 0.05 M TRIS buffer, pH 8.0.
Preservative	None
Storage	The antibody should be stored at 2-8°C. For storage longer than one year, the solution may be frozen at -20°C. Repeated freezing and thawing is not recommended. The solution may be frozen in aliquots if necessary. The antibody has a shelf-life of 2 years after date of purchase.

BACKGROUND

Introduction	PFBV is mainly occurring in Europe. It is restricted to geraniums (<i>Pelargonium zonale</i> and <i>Pelargonium peltatum</i>). It is transmitted by inoculation of sap. It causes flower-breaking in some pelargonium seedling lines, but is symptomless in many cultivars.
Keywords	Pelargonium flower-break virus; PFBV