



# Anti-Vibrio Cholera Toxin A+B polyclonal antibody (CABT-BL8725)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Target</b>	Cholera Toxin
<b>Isotype</b>	IgG
<b>Source/Host</b>	Goat
<b>Species Reactivity</b>	Cholera
<b>Purification</b>	Antigen Affinity Purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, IP, WB
<b>Format</b>	Liquid
<b>Size</b>	0.1 mg
<b>Buffer</b>	Phosphate buffered saline, pH 7.2
<b>Preservative</b>	None
<b>Storage</b>	4°C for 2 years