



Anti-CD40LG monoclonal antibody, clone TB-EB5 [FITC] (CABT-BL947)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	TB-EB5
Conjugate	FITC
Applications	FC
Preservative	None
Storage	-20°C. Avoid freeze/thaw

BACKGROUND

Introduction The protein encoded by this gene is expressed on the surface of T cells. It regulates B cell function by engaging CD40 on the B cell surface. A defect in this gene results in an inability to undergo immunoglobulin class switch and is associated with hyper-IgM syndrome. [provided by RefSeq, Jul 2008]

GENE INFORMATION

Entrez Gene ID	959
Protein Refseq	NP_000065
UniProt ID	P29965