



Mouse Anti-Methyl Mercury Chloride monoclonal antibody, clone CH3HgCl (CABT-L2514)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Methyl Mercury Chloride
Immunogen	Methyl Mercury Chloride with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	CH3HgCl
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	-20°C. Avoid Freeze/Thaw Cycles
Warnings	For research use only and not for use in diagnostic or therapeutic procedures.

BACKGROUND

Introduction Methylmercury (sometimes methyl mercury) is an organometallic cation with the formula [CH₃Hg]⁺. It is the major source of organic mercury for all humans.