



# Mouse Anti-Maduramicin monoclonal antibody, clone MAD (CABT-L2535)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	Maduramicin
<b>Immunogen</b>	Maduramicin with carrier protein
<b>Isotype</b>	IgG
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	N/A
<b>Clone</b>	MAD
<b>Purification</b>	Protein A purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, LFIA
<b>Format</b>	Liquid
<b>Concentration</b>	Lot specific
<b>Size</b>	1 mg
<b>Buffer</b>	PBS
<b>Preservative</b>	None
<b>Storage</b>	-20°C. Avoid Freeze/Thaw Cycles
<b>Warnings</b>	For research use only and not for use in diagnostic or therapeutic procedures.

## BACKGROUND

<b>Introduction</b>	Maduramicin is a monoglycoside polyether derived from the fungus <i>Actinomadura rubra</i> .
<b>Keywords</b>	Maduramicin;Maduramycin;Maduramicin Ammonium;Maduramycin Ammonium;Cygro;